

<b>LYNX EVO</b>									
	<b>Zoom Ratio – 10x</b>			<b>Zoom Ratio – 15x</b>			<b>Zoom Ratio – 20x</b>		
				<b>1.5 Multiplier</b>			<b>2.0 Multiplier</b>		
<b>Objective Lenses</b>	<b>Mag X</b>	<b>WD (mm)</b>	<b>FOV (mm)</b>	<b>Mag X</b>	<b>WD (mm)</b>	<b>FOV (mm)</b>	<b>Mag X</b>	<b>WD (mm)</b>	<b>FOV (mm)</b>
<b>0.45x</b>	2.7x – 27x	176	55	2.7x – 40x	176	55	2.7x – 54x	176	55
<b>0.62x</b>	3.7x – 37x	128	40	3.7x – 55x	128	40	3.7x – 74x	128	40
<b>1.0x</b>	6x – 60x	75	25	6x – 90x	75	25	6x – 120x	75	25
<b>1.5x</b>	9x – 90x	42	16	9x – 135x	42	16	9x – 180x	42	16
<b>2.0x</b>	12x – 120x	29	12	12x – 180x	29	12	12x – 240x	29	12

<b>LYNX EVO ~ Lite</b>									
	<b>Direct</b>			<b>1.5 Multiplier</b>			<b>2.0 Multiplier</b>		
<b>Objective Lenses</b>	<b>Mag X</b>	<b>WD (mm)</b>	<b>FOV (mm)</b>	<b>Mag X</b>	<b>WD (mm)</b>	<b>FOV (mm)</b>	<b>Mag X</b>	<b>WD (mm)</b>	<b>FOV (mm)</b>
<b>0.45x</b>	8.5x	175	17.5	11x	175	11.6	17x	175	8.9
<b>0.5x</b>	10x	159	15.3	15x	159	10.2	20x	159	7.6
<b>0.62x</b>	12x	125	12.9	18x	125	8.6	24x	125	6.4
<b>1.0x</b>	20x	75	7.6	30x	75	5.1	40x	75	3.8
<b>1.5x</b>	30x	42	5.1	45x	42	3.4	60x	42	2.6
<b>2.0x</b>	40x	29	3.8	60x	29	2.55	80x	29	1.9